

Sieve Diameters and Frame Materials

Diameter	Height	Depth to Mesh or Plate	Frame Material
38	Full	19mm	Br or SS
100	Full	40mm	Br or SS
100	Half	20mm	Br or SS
150	Full	38mm	SS
200	Full	50mm	Br or SS
200	Half	25mm	Br or SS
250	Full	60mm	SS
300	Full	75mm	Br,SS or PS
300	Half	40mm	Br,SS or PS
315	Full	75mm	SS
350	Full	60mm	Br or SS
400	Full	65mm	Br,SS or PS
450	Full	100mm	SS or PS
3 in	Full	1¼ in	Br or SS
8 in	Full	2 in	Br or SS
8 in	Half	1 in	Br or SS
12 in	Full	3 in	Br,SS or PS
12 in	Half	1½ in	Br,SS or PS
18 in	Full	3½ in	SS or PS